

Product datasheet for **GA113242**

Human Junctional Adhesion Molecule C (JAM3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	JAM3 CRISPRa kit - CRISPR gene activation of human junctional adhesion molecule 3
Symbol:	JAM3
UniProt ID:	Q9BX67
Synonyms:	JAM-2; JAM-3; JAM-C; JAMC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)