

Product datasheet for **GA113249**

Human TLCD3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLCD3B CRISPRa kit - CRISPR gene activation of human TLC domain containing 3B
Symbol:	TLCD3B
Components:	GA113249G1, TLCD3B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113249G2, TLCD3B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113249G3, TLCD3B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q71RH2</u>
Synonyms:	FP1188
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)