

Product datasheet for **GA113255**

Human Histone H2A Bbd (H2AFB3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	H2AFB3 CRISPRa kit - CRISPR gene activation of human H2A histone family member B3
Symbol:	H2AB3
UniProt ID:	<u>P0C5Z0</u>
Synonyms:	H2ABBD; H2AFB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)