

Product datasheet for **GA113268**

Human **KCNK16** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK16 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 16
Symbol:	KCNK16
UniProt ID:	Q96T55
Synonyms:	K2p16.1; TALK-1; TALK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)