

Product datasheet for **GA113279**

Human ATAD3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATAD3B CRISPRa kit - CRISPR gene activation of human ATPase family AAA domain containing 3B
Symbol:	ATAD3B
Components:	GA113279G1, ATAD3B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113279G2, ATAD3B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113279G3, ATAD3B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5T9A4
Synonyms:	AAA-TOB3; TOB3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)