


Product datasheet for **GA113280**

Human TAF3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAF3 CRISPRa kit - CRISPR gene activation of human TATA-box binding protein associated factor 3
Symbol:	TAF3
Components:	GA113280G1, TAF3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113280G2, TAF3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113280G3, TAF3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5VWG9</u>
Synonyms:	TAF140; TAFII-140; TAFII140
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)