

## Product datasheet for **GA113286**

### Human TBC1D10A activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | TBC1D10A CRISPRa kit - CRISPR gene activation of human TBC1 domain family member 10A                                                                                                                                                                                                                                |
| Symbol:       | TBC1D10A                                                                                                                                                                                                                                                                                                            |
| Components:   | GA113286G1, TBC1D10A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113286G2, TBC1D10A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113286G3, TBC1D10A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9BXI6</a>                                                                                                                                                                                                                                                                                              |
| Synonyms:     | dj130H16.1; dj130H16.2; EPI64; TBC1D10                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)