

Product datasheet for **GA113313**

Human **STK40** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	STK40 CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 40
Symbol:	STK40
UniProt ID:	Q8N2I9
Synonyms:	SgK495; SHIK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)