

Product datasheet for **GA113317**

Human LRMDA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRMDA CRISPRa kit - CRISPR gene activation of human leucine rich melanocyte differentiation associated
Symbol:	LRMDA
Components:	GA113317G1, LRMDA gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113317G2, LRMDA gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113317G3, LRMDA gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H2I8
Synonyms:	CDA017
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

