

Product datasheet for **GA113350**

Human C1orf49 (TEX35) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEX35 CRISPRa kit - CRISPR gene activation of human testis expressed 35
Symbol:	TEX35
UniProt ID:	<u>Q5T0J7</u>
Synonyms:	C1orf49; TSC24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)