

Product datasheet for **GA113353**

Human PLEKHN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHN1 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing N1
Symbol:	PLEKHN1
UniProt ID:	Q494U1
Synonyms:	RP11-5407.7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)