

Product datasheet for **GA113366**

Human ZRANB3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZRANB3 CRISPRa kit - CRISPR gene activation of human zinc finger RANBP2-type containing 3
Symbol:	ZRANB3
Components:	GA113366G1, ZRANB3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113366G2, ZRANB3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113366G3, ZRANB3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5FWF4</u>
Synonyms:	4933425L19Rik; AH2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)