

Product datasheet for **GA113372**

Human C4orf17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C4orf17 CRISPRa kit - CRISPR gene activation of human chromosome 4 open reading frame 17
Symbol:	C4orf17
Components:	GA113372G1, C4orf17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113372G2, C4orf17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113372G3, C4orf17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q53FE4
Synonyms:	DKFZp434G072
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)