

Product datasheet for **GA113378**

Human ZNF394 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF394 CRISPRa kit - CRISPR gene activation of human zinc finger protein 394
Symbol:	ZNF394
UniProt ID:	Q53GI3
Synonyms:	ZKSCAN14; ZSCAN46
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)