

## Product datasheet for **GA113380**

### Human ATRIP activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | ATRIP CRISPRa kit - CRISPR gene activation of human ATR interacting protein |
| Symbol:       | ATRIP   |
| UniProt ID:   | <u><a href="#">Q8WXE1</a></u>   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)     |

#### Product images:



[View online >](#)