

Product datasheet for **GA113383**

Human CEP78 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CEP78 CRISPRa kit - CRISPR gene activation of human centrosomal protein 78
Symbol:	CEP78
Components:	GA113383G1, CEP78 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113383G2, CEP78 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113383G3, CEP78 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5JTW2</u>
Synonyms:	C9orf81; IP63
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)