

Product datasheet for **GA113393**

Human ZNF644 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF644 CRISPRa kit - CRISPR gene activation of human zinc finger protein 644
Symbol:	ZNF644
Components:	GA113393G1, ZNF644 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113393G2, ZNF644 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113393G3, ZNF644 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H582
Synonyms:	BM-005; MYP21; NatF; ZEP-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)