

Product datasheet for **GA113403**

Human C19orf44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C19orf44 CRISPRa kit - CRISPR gene activation of human chromosome 19 open reading frame 44
Symbol:	C19orf44
Components:	GA113403G1, C19orf44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113403G2, C19orf44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113403G3, C19orf44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6X5</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)