


Product datasheet for **GA113409**

Human SLC49A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC49A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 49 member 3
Symbol:	SLC49A3
Components:	GA113409G1, SLC49A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113409G2, SLC49A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113409G3, SLC49A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UXD7</u>
Synonyms:	LP2561
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)