


Product datasheet for **GA113413**

Human TMEM164 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM164 CRISPRa kit - CRISPR gene activation of human transmembrane protein 164
Symbol:	TMEM164
Components:	GA113413G1, TMEM164 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113413G2, TMEM164 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113413G3, TMEM164 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5U3C3</u>
Synonyms:	bB360B22.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)