

Product datasheet for **GA113418**

Human SETD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SETD3 CRISPRa kit - CRISPR gene activation of human SET domain containing 3, actin histidine methyltransferase
Symbol:	SETD3
Components:	GA113418G1, SETD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113418G2, SETD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113418G3, SETD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86TU7</u>
Synonyms:	C14orf154
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)