

Product datasheet for **GA113423**

Human ANKRD20A1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ANKRD20A1 CRISPRa kit – CRISPR gene activation of human ankyrin repeat domain 20 family member A1
Symbol:	ANKRD20A1
Components:	GA113423G1, ANKRD20A1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113423G2, ANKRD20A1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113423G3, ANKRD20A1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TYW2</u>
Synonyms:	ANKRD20A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

