


Product datasheet for **GA113424**

Human ZNF541 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ZNF541 CRISPRa kit - CRISPR gene activation of human zinc finger protein 541 |
| Symbol: | ZNF541 |
| UniProt ID: | Q9H0D2 |
| Synonyms: | DKFZp434I1930; MGC119518; MGC119519; MGC119521 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)