


## Product datasheet for **GA113428**

### Human WDR24 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | WDR24 CRISPRa kit - CRISPR gene activation of human WD repeat domain 24                                                                                                                                                                                                                                    |
| Symbol:       | WDR24                                                                                                                                                                                                                                                                                                      |
| Components:   | GA113428G1, WDR24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113428G2, WDR24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113428G3, WDR24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q96S15</a>                                                                                                                                                                                                                                                                                     |
| Synonyms:     | C16orf21; JFP7                                                                                                                                                                                                                                                                                             |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)