

Product datasheet for **GA113451**

Human GARNL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GARNL3 CRISPRa kit - CRISPR gene activation of human GTPase activating Rap/RanGAP domain like 3
Symbol:	GARNL3
Components:	GA113451G1, GARNL3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113451G2, GARNL3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113451G3, GARNL3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5VWW2
Synonyms:	bA356B19.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)