


Product datasheet for **GA113455**

Human Synaptotagmin 3 (SYT3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SYT3 CRISPRa kit - CRISPR gene activation of human synaptotagmin 3
Symbol:	SYT3
UniProt ID:	<u>Q9BQG1</u>
Synonyms:	SytIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)