

Product datasheet for **GA113472**

Human **SLC25A33** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A33 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 33
Symbol:	SLC25A33
UniProt ID:	Q9BSK2
Synonyms:	BMSC-MCP; PNC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)