

## Product datasheet for **GA113483**

### Human Transmembrane Protein 175 (TMEM175) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM175 CRISPRa kit - CRISPR gene activation of human transmembrane protein 175
Symbol:	TMEM175
UniProt ID:	<a href="#">Q9BSA9</a>
Synonyms:	hTMEM175
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)