


Product datasheet for **GA113488**

Human Morg1 (WDR83) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR83 CRISPRa kit - CRISPR gene activation of human WD repeat domain 83
Symbol:	WDR83
Components:	GA113488G1, Morg1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113488G2, Morg1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113488G3, Morg1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BRX9
Synonyms:	MORG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)