

Product datasheet for **GA113502**

Human ZNF397 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF397 CRISPRa kit - CRISPR gene activation of human zinc finger protein 397
Symbol:	ZNF397
Components:	GA113502G1, ZNF397 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113502G2, ZNF397 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113502G3, ZNF397 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8NF99</u>
Synonyms:	ZNF47; ZSCAN15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)