

Product datasheet for **GA113508**

Human TMEM107 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM107 CRISPRa kit - CRISPR gene activation of human transmembrane protein 107
Symbol:	TMEM107
Components:	GA113508G1, TMEM107 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113508G2, TMEM107 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113508G3, TMEM107 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6UX40</u>
Synonyms:	GRVS638; PRO1268
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)