


Product datasheet for **GA113510**

Human LSMD1 (NAA38) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NAA38 CRISPRa kit - CRISPR gene activation of human N(alpha)-acetyltransferase 38, NatC auxiliary subunit
Symbol:	NAA38
Components:	GA113510G1, LSMD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113510G2, LSMD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113510G3, LSMD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BRA0</u>
Synonyms:	LSMD1; PFAAP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)