

## Product datasheet for **GA113522**

### Human ZNF414 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | ZNF414 CRISPRa kit – CRISPR gene activation of human zinc finger protein 414                                                                                                                                                                                                                                  |
| <b>Symbol:</b>       | ZNF414                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA113522G1, ZNF414 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113522G2, ZNF414 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA113522G3, ZNF414 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <a href="#">Q96IQ9</a>                                                                                                                                                                                                                                                                                        |
| <b>Synonyms:</b>     | ZFP414                                                                                                                                                                                                                                                                                                        |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

