

Product datasheet for **GA113528**

Human Transmembrane Protein 101 (TMEM101) activation kit by CRISPRa


Product data:

Product Type:	CRISPRa Kits
Description:	TMEM101 CRISPRa kit - CRISPR gene activation of human transmembrane protein 101
Symbol:	TMEM101
Components:	GA113528G1, Transmembrane Protein 101 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113528G2, Transmembrane Protein 101 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113528G3, Transmembrane Protein 101 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077

UniProt ID: [Q96IK0](#)

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)