

Product datasheet for GA113532

Human HPS3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HPS3 CRISPRa kit - CRISPR gene activation of human HPS3 biogenesis of lysosomal organelles complex 2 subunit 1
Symbol:	HPS3
UniProt ID:	Q969F9
Synonyms:	BLOC2S1; SUTAL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)