


## Product datasheet for **GA113549**

### Human DOT1L activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	DOT1L CRISPRa kit - CRISPR gene activation of human DOT1 like histone lysine methyltransferase
Symbol:	DOT1L
UniProt ID:	<u><a href="#">Q8TEK3</a></u>
Synonyms:	DOT1; KMT4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)