


## Product datasheet for **GA113592**

### Human PGBD1 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | PGBD1 CRISPRa kit - CRISPR gene activation of human piggyBac transposable element derived 1 |
| Symbol:       | PGBD1   |
| UniProt ID:   | <a href="#">Q96JS3</a>  |
| Synonyms:     | dJ874C20.4; HUCEP-4; SCAND4   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                     |

#### Product images:



[View online >](#)