

Product datasheet for **GA113604**

Human ZBTB37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZBTB37 CRISPRa kit - CRISPR gene activation of human zinc finger and BTB domain containing 37
Symbol:	ZBTB37
Components:	GA113604G1, ZBTB37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113604G2, ZBTB37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113604G3, ZBTB37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TC79</u>
Synonyms:	D430004I08Rik; ZNF908
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)