


Product datasheet for **GA113616**

Human TNRC18 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNRC18 CRISPRa kit - CRISPR gene activation of human trinucleotide repeat containing 18
Symbol:	TNRC18
Components:	GA113616G1, TNRC18 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113616G2, TNRC18 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113616G3, TNRC18 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q15417</u>
Synonyms:	CAGL79; TNRC18A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)