

Product datasheet for **GA113622**

Human **IL1F10** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	IL1F10 CRISPRa kit - CRISPR gene activation of human interleukin 1 family member 10
Symbol:	IL1F10
Components:	GA113622G1, IL1F10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113622G2, IL1F10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113622G3, IL1F10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WWZ1</u>
Synonyms:	FIL1-theta; FKSG75; IL-1HY2; IL-38; IL1-theta; IL1HY2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)