

Product datasheet for **GA113628**

Human LCE3D activation kit by CRISPRa

Product data:

| | |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | LCE3D CRISPRa kit - CRISPR gene activation of human late cornified envelope 3D |
| Symbol: | LCE3D |
| Components: | GA113628G1, LCE3D gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113628G2, LCE3D gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113628G3, LCE3D gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>Q9BYE3</u> |
| Synonyms: | LEP16; SPRL6A; SPRL6B |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)