

Product datasheet for **GA113633**

Human GLYR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GLYR1 CRISPRa kit - CRISPR gene activation of human glyoxylate reductase 1 homolog
Symbol:	GLYR1
Components:	GA113633G1, GLYR1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113633G2, GLYR1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113633G3, GLYR1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q49A26</u>
Synonyms:	BM045; HIBDL; N-PAC; NP60
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)