

Product datasheet for **GA113674**

Human PSRC1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSRC1 CRISPRa kit - CRISPR gene activation of human proline and serine rich coiled-coil 1
Symbol:	PSRC1
UniProt ID:	Q6PGN9
Synonyms:	DDA3; FP3214
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)