

Product datasheet for **GA113675**

Human FAPP2 (PLEKHA8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA8 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A8
Symbol:	PLEKHA8
Components:	GA113675G1, FAPP2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113675G2, FAPP2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113675G3, FAPP2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96JA3
Synonyms:	FAPP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)