


Product datasheet for **GA113689**

Human FKSG24 (MPV17L2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MPV17L2 CRISPRa kit - CRISPR gene activation of human MPV17 mitochondrial inner membrane protein like 2
Symbol:	MPV17L2
Components:	GA113689G1, FKSG24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113689G2, FKSG24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113689G3, FKSG24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q567V2</u>
Synonyms:	FKSG24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)