

## Product datasheet for **GA113703**

### Human PPAPDC3 (PLPP7) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PLPP7 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase 7 (inactive)
Symbol:	PLPP7
UniProt ID:	<a href="#">Q8NBV4</a>
Synonyms:	C9orf67; NET39; PPAPDC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)