

Product datasheet for **GA113709**

Human SFT2C (SFT2D3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SFT2D3 CRISPRa kit - CRISPR gene activation of human SFT2 domain containing 3
Symbol:	SFT2D3
Components:	GA113709G1, SFT2C gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113709G2, SFT2C gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113709G3, SFT2C gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q587I9</u>
Synonyms:	MGC5391
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)