

Product datasheet for **GA113711**

Human ATP5MK activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATP5MD CRISPRa kit - CRISPR gene activation of human ATP synthase membrane subunit DAPIT
Symbol:	ATP5MK
Components:	GA113711G1, ATP5MK gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113711G2, ATP5MK gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113711G3, ATP5MK gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96IX5
Synonyms:	bA792D24.4; DAPIT; HCVFTP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)