

Product datasheet for **GA113713**

Human ZNF496 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF496 CRISPRa kit - CRISPR gene activation of human zinc finger protein 496
Symbol:	ZNF496
Components:	GA113713G1, ZNF496 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113713G2, ZNF496 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113713G3, ZNF496 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96IT1</u>
Synonyms:	NIZP1; ZFP496; ZKSCAN17; ZSCAN49
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)