

## Product datasheet for **GA113736**

### Human **RPUSD4** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	RPUSD4 CRISPRa kit - CRISPR gene activation of human RNA pseudouridine synthase D4
Symbol:	RPUSD4
Components:	GA113736G1, RPUSD4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA113736G2, RPUSD4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA113736G3, RPUSD4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q96CM3</a></u>
Synonyms:	FLJ14494
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)